DIL-Series

DIL1-B1

404450305887

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Intended use

The DIL series instruments are used for the endoscopic dilation of strictures in the gastrointestinal tract. The inflation pump fills the dilation balloon for dilation.

Product features

incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.


ATTENTION! Use only sterile water, saline solution or contrast agent solutions as the inflation medium. Do not use air or gases for inflation purposes, as these reduce the effectiveness of the balloon.

Instructions for use

Insert the connecting tube into the container filled with fluid, press the yellow shaft and pull back the black handle. When the plunger is sufficiently filled with liquid, release the yellow shaft. Then purge the system of air. Attach the connecting tube to the three-way tap on the dilation balloon. Press the yellow shaft and press the fluid into the balloon using the black handle. Release the yellow shaft. To build up pressure, turn the black handle clockwise. To reduce pressure, turn the black handle anti-clockwise, until you have reached 2-3 bar. Now press the yellow shaft and pull the black handle back until the balloon is fully deflated. Remove the dilation instrument on the three-way tap from the dilation balloon.

After completion of the surgical procedure

Instruments intended for single use, including packaging, must be disposed of in accordance with the relevant applicable hospital and administrative guidelines and in accordance with the relevant applicable legal regulations.

Qualification of user

Use of these instruments requires comprehensive knowledge of the technical principles, clinical applications and risks of gastrointestinal endoscopy. The instruments should only be used by or under the supervision of physicians who are adequately trained and experienced in endoscopic techniques.

General instructions

Use this instrument only and exclusively for the purposes described in these instructions.

ATTENTION! Instruments marked with this symbol are intended for single use only and are sterilised with ethylene oxide.

Sterile instruments can be used immediately. Before use, check the "Use by" date on the packaging, as sterile instruments are only allowed to be used up to this date.

Do NOT use the instrument if the sterile packaging has tears or perforations, sealing is not ensured or moisture has penetrated. All medwork instruments should be stored in a dry place that is protected from light. Keep all instructions for use in a safe and easily accessible place.

medwork instruments that are labelled as intended for single use only must not be reprocessed, resterilised or reused. Reuse, reprocessing or resterilisation can alter product properties, resulting in functional failure that can endanger the patient's health or lead to disease, injury or death. Reuse, reprocessing or resterilisation furthermore entails the risk of contamination of the patient or the instrument, as well as the risk of cross contamination, including the transmission of infectious diseases. Contamination of the instrument can result in illness, injury or death of the patient.

Indications

Strictures in the upper and lower gastrointestinal tract, gastroenteric anastomotic strictures, benign biliary strictures, achalasia, Crohn's strictures, biliary and pancreatic stones.

Contraindications

The contraindications for the DIL series correspond to the specific contraindications for oesophagogastroduodenoscopy and endoscopic retrograde cholangiopancreatography.

Possible complications

Allergic reaction to contrast agents, bleeding, perforation, sepsis/infection, pancreatitis, cholangitis.

Precautions

In order to ensure that the examination proceeds smoothly, the diameter of the endoscope working channel and the diameter of the instruments must be matched to one another.

After removal from the packaging, check the instruments for correct function, kinks, fractures, rough surfaces, sharp edges and protrusions. If you discover any damage or malfunction on the instruments, do NOT use them and inform your responsible field service contact or our branch office.

Any person who markets or uses medical devices is required to report any serious

List of all icons used on medwork products

	Date of manufacture		Use by
	Follow the instructions for use		Caution near pacemakers
	Applied part Type BF		Do not use if packaging is damaged
	Batch code		Not suitable for lithotripsy
	Item number		Contains natural rubber latex
	Packaging unit		Gastroscopy
	Do not reuse		Colonoscopy
	Sterilised using ethylene oxide		Enteroscopy
	Non-sterile		ERCP
	Manufacturer		Medical device
	Sterile barrier system		

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